FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Material Safety Data Sheet (MSDS)

Section 1 — Chemical Product and Company Identification

o-Nitrobenzaldehyde

Flinn Scientific, Inc. P.O. Box 219 Batavia, IL 60510 (800) 452-1261 CHEMTREC Emergency Phone Number: (800) 424-9300

Section 2 - Composition, Information on Ingredients

ortho-Nitrobenzaldehyde

CAS#: 552-89-6

Section 3 — Hazards Identification

Light yellow or buff crystalline powder. Moderately toxic by ingestion. Eye irritant. Combustible solid.

Health-2 Flammability-1 Reactivity-1 Exposure-1 Storage-0

FLINN AT-A-GLANCE

0 is low hazard, 3 is high hazard

Section 4 — First Aid Measures

Call a physician, seek medical attention for further treatment, observation and support after first aid. Inhalation: Remove to fresh air at once. If breathing has stopped give artificial respiration immediately. Eye: Immediately flush with fresh water for 15 minutes. External: Wash continuously with fresh water for 15 minutes.

Internal: Rinse out mouth, give 1 to 2 cups of water or milk, induce vomiting. Call a physician or poison control at once.

Section 5 — Fire Fighting Measures

Combustible solid. When heated to decomposition, may emit toxic fumes of NOx. Fire Fighting Instructions: Use triclass, dry chemical fire extinguisher. Firefighters should wear PPE and SCBA with full facepiece operated in positive pressure mode.

Section 6 — Accidental Release Measures

Restrict unprotected personnel from area. Remove all ignition sources and water. Sweep up, place in sealed bag or container and dispose. Ventilate area and wash spill site after material pickup is complete. See Sections 8 and 13 for further information.

Section 7 — Handling and Storage

Flinn Suggested Chemical Storage Pattern: Organic #3. Store with hydrocarbons, esters and aldehydes. Store in a cool, dry place.

Section 8 — Exposure Controls, Personal Protection

Avoid contact with eyes, skin and clothing. Wear chemical splash goggles, chemical-resistant gloves and chemical-resistant apron. Use exhaust ventilation to keep airborne concentrations low.

PAGE 1 OF 2

MSDS #: 553.30 Revision Date: April 2, 2001

NFPA CODE None established

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

o-Nitrobenzaldehyde

Material Safety Data Sheet (MSDS)

MSDS #: 553.30 Revision Date: April 2, 2001

Section 9 — Physical and Chemical Properties

Light yellow or buff powder. Solubility: Slightly soluble in water; soluble in alcohol, ether, and acetone. Formula: C7H5O3N Formula Weight: 151.12 Melting Point: 42-44 C

Section 10 — Stability and Reactivity

Avoid contact with strong oxidizers and bases. Shelf life: Indefinite.

Section 11 — Toxicological Information

Acute effects: Irritant Chronic effects: N.A. Target organs: N.A. ORL-MUS LD50: 600 mg/kg IHL-RAT LC50: N.A. SKN-RBT LD50: N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

Section 12 — Ecological Information

Data not yet available.

Section 13 — Disposal Considerations

Please consult with state and local regulations before disposal. Flinn Suggested Disposal Method #2 is one option.

Section 14 — Transport Information

Shipping Name: Not regulated Hazard Class: N/A UN Number: N/A

N/A = Not applicable

Section 15 — Regulatory Information

TSCA-listed, EINECS-listed (209-025-3).

Section 16 — Other Information

Consult your copy of the Flinn Scientific Catalog/Reference Manual for additional information about laboratory chemicals. This Material Safety Data Sheet (MSDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. Flinn Scientific Inc. assumes no legal responsibility for use or reliance upon this data.

FLINN SCIENTIFIC INC.

"Your Safer Source for Science Supplies"

Questions on Chemical Disposal or Storage?--Call Flinn

> flinn@flinnsci.com www.flinnsci.com P.O. Box 219 Batavia IL 60510 (800) 452-1261 Fax (866) 452-1436